



QUALITY CONTROL RANGES

Mouse MMP Magnetic Bead Panel 1
Kit Catalog No. MMMP1MAG-79K

Mouse MMP Magnetic Bead Panel 2
Kit Catalog No. MMMP2MAG-79K

Mouse MMP Magnetic Bead Panel 3
Kit Catalog No. MMMP3MAG-79K

Control Catalog No. MMMP-6079
Lot No. MMMP-111 and MMMP-211

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
MMMP2	Control 1	574 - 1193	pg/mL
	Control 2	5996 - 12453	pg/mL
MMMP3	Control 1	32 - 66	pg/mL
	Control 2	329 - 683	pg/mL
MMMP8	Control 1	148 - 307	pg/mL
	Control 2	1529 - 3175	pg/mL
MMMP9	Control 1	613 - 1273	pg/mL
	Control 2	6525 - 13551	pg/mL
MMMP12	Control 1	56 - 117	pg/mL
	Control 2	631 - 1311	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format.
MMMP-6079: MMMP-111.211 (04 November 2025) • [SigmaAldrich.com/support](https://www.sigmaldrich.com/support) •
[SigmaAldrich.com](https://www.sigmaldrich.com)