


Anti-ASAP1/Centaurin B4 (ED-20) (A 4227)



Immunoblot

Mouse brain (Lanes A, C) and rat brain (Lane B) cytosolic fractions were separated on SDS-PAGE and blotted with Rabbit Anti-ASAP1 (ED-20) (Product No. A 4227). The Antibody was developed with Goat Anti-Rabbit IgG Peroxidase (Product No. A 0545) and chemiluminescent substrate.

Lanes A, B: Antibody dilution 1: 2,500

Lane C: Antibody dilution 1: 2,500 + 10 μ g/ml ASAP1 immunizing peptide.