

M0553 Middlebrook ADC Growth Supplement

An enrichment supplement recommended for cultivation of Mycobacteria.

Composition

Bovine albumin fraction V 2.5 g

Dextrose 1 g

Catalase 0.0015 g

Distilled water 50 ml

Store at 4°C.

Directions

Warm up the refrigerated supplement to 45-50°C. Aseptically add the contents of the vial to 450 ml of sterile Middlebrook 7H9 Broth Base (Cat. No. M0178). Mix well and pour into sterile screw capped tubes.

