besign for sustainability (bfs) scorecard

Results

With our DfS scorecard, we drive sustainability improvement during the product development process through multiple product sustainability criteria divided into seven impact areas.



Milli-Q[®] IX 7010 Pure Water Purification System

Ultimate pure water solution - less electricity, more compact and mercury-free UV lamps



Im	pa	ct	aı	e	as
	?				

MATERIALS

15% system size reduction but 4% weight increase due to more efficient but heavier power supply 60% reduction of plastic use in purification cartridges thanks to extended lifespan and more compact design



SUPPLIERS & MANUFACTURING

No change compared to baseline product in consideration of our DfS criteria



PACKAGING

For purification cartridges' packaging, 32% weight and 30% volume reduction compared to initial design due to new protective insert design and twin box concept. System's corrugated box with sustainable forestry certification. Protective inserts made of 100% bio-based polyethylene (PE) foam from sugar cane waste for the system and made of PE foam with a minimum of 50% recycled content for the E-POD® water remote dispensers



ENERGY & EMISSIONS

32% reduction of system electricity consumption thanks to optimization of components and processes as well as an embedded improved "Lab Close" mode



WATER

7% water reduction during system use thanks to new periodic EDI rinsing
Reduction of water consumption during system quality control thanks to the implementation of a new impermeability
test by air instead of by water



USABILITY & INNOVATION

Innovative bactericidal lamp with mercury-free UVC LED technology Easy and quick cartridge replacement thanks to unique patented twist and lock technology



CIRCULAR ECONOMY

The system is designed to be easily maintainable and have a lower replacement frequency of purification consumables. Implementation of a drain cap enables to purge about 44% of the water from the purification cartridges, representing about 25% of the cartridges' total weight prior to disposal.

Baseline product: ZRXV010T0 - Elix® Advantage 10