



QUALITY CONTROL RANGES

Human Pancreatic Polypeptide ELISA Kit Catalog #EZHPP40K

Human PP Standard, Lyophilized Catalog # E8040K

STANDARD	HUMAN PP ELISA VALUE
Lot Number: EHPP-13K	2818 pg/mL

The exact value of lyophilized standard is shown above.
Reconstitute with **0.5mL** distilled or deionized water.

Quality Controls 1 and 2, Lyophilized Control Catalog # E6040K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HPP-113	119 - 247	pg/mL
QC 2	HPP-213	595 - 1235	pg/mL

Reconstitute each vial with **0.5mL** distilled or deionized water.

E6040K: HPP-113.213 (02 April 2024)
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)
www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA