

# Product Information

## Pluronic® F68 Solution 10%

100.0 g/L Pluronic® F68

CATALOG NO. 59915C

### Description

In low-serum or serum-free conditions, cells in suspension are exposed to destructive mechanical stresses. The addition of 0.1% Pluronic® F68 to culture media protects the cells by reducing shear. This 10% solution should be diluted 1:100 prior to use. Concentrations higher than 0.1% can be cytotoxic to some cell types.

### Precautions

Use aseptic technique when handling or using this supplement. This product is for further manufacturing use. THIS PRODUCT IS NOT INTENDED FOR HUMAN OR THERAPEUTIC USE.

### Storage

To effectively preserve the integrity of Pluronic® F68, it should be stored at the recommended temperature and protected from light. For stability and optimal performance, this product should be stored at 15 to 30 C. Do not use after the expiration date.

### Indications of Deterioration

Supplement should be clear of particulates and flocculent material. Do not use if supplement is cloudy or contains precipitate. Other evidence of deterioration may include color change as well as degradation of physical or performance characteristics.

### Characteristics

#### Appearance

Clear colorless solution

#### Endotoxin (LAL)

≤ 1.0 EU/mL

#### Osmolality (as supplied)

73 - 113 mOsm/kg H<sub>2</sub>O

#### pH (as supplied)

6.5 - 7.0

#### Sterility

No microbial growth detected

#### United States

SAFC Biosciences, Inc.  
13804 W. 107th Street  
Lenexa, Kansas 66215  
USA  
Phone +1 913-469-5580  
Toll free-USA 1 800-255-6032  
Fax +1 913-469-5584  
E-mail info-na@sial.com

#### Europe

SAFC Biosciences Ltd.  
Smeaton Road, West Portway  
Andover, Hampshire SP10 3LF  
UNITED KINGDOM  
Phone +44 (0)1264-333311  
Fax +44 (0)1264-332412  
E-mail info-eu@sial.com

#### Asia Pacific

SAFC Biosciences Pty. Ltd.  
18-20 Export Drive  
Brooklyn, Victoria 3025  
AUSTRALIA  
Phone +61 (0)3-9362-4500  
Toll free-AUS 1 800-200-404  
Fax +61 (0)3-9315-1656  
E-mail info-ap@sial.com

#### Warranty, Limitation of Remedies

SAFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC Biosciences expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

Pluronic® is a registered trademark of BASF Corporation.

© 2006 SAFC Biosciences, Inc.

Issued December 2006 P59915  
0103 0905 1005 0406 0906

#### United States

SAFC Biosciences, Inc.  
13804 W. 107th Street  
Lenexa, Kansas 66215  
USA  
Phone +1 913-469-5580  
Toll free-USA 1 800-255-6032  
Fax +1 913-469-5584  
E-mail info-na@sial.com

#### Europe

SAFC Biosciences Ltd.  
Smeaton Road, West Portway  
Andover, Hampshire SP10 3LF  
UNITED KINGDOM  
Phone +44 (0)1264-333311  
Fax +44 (0)1264-332412  
E-mail info-eu@sial.com

#### Asia Pacific

SAFC Biosciences Pty. Ltd.  
18-20 Export Drive  
Brooklyn, Victoria 3025  
AUSTRALIA  
Phone +61 (0)3-9362-4500  
Toll free-AUS 1 800-200-404  
Fax +61 (0)3-9315-1656  
E-mail info-ap@sial.com