



# QUALITY CONTROL RANGES

## CANINE LEPTIN ELISA KIT CATALOG # EZCL-31K

### CANINE LEPTIN STANDARD, LYOPHILIZED CATALOG # E8031-K

STANDARD	CANINE LEPTIN VALUE
Lot Number: 831-12K	54 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6031-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	CL-112	1.8 – 3.7	ng/mL
QC 2	CL-212	11 – 23	ng/mL

E6031-K: CL-112.212 (15 November 2023)  
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

[www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.  
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA