



Quality Control Ranges
MILLIPLEX® Human Exosome Characterization Magnetic Bead Panel
Catalog # HEXSM-170K
Control Catalog # HEXSM-6170 Lot # HEXSM-101.201

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
CALR	Control 1	1381 - 2868	pg/mL
	Control 2	5519 - 11462	pg/mL
CD9	Control 1	24 - 50	ng/mL
	Control 2	116 - 241	ng/mL
CD63	Control 1	6862 - 14253	pg/mL
	Control 2	27052 - 56185	pg/mL
CD81	Control 1	28 - 57	ng/mL
	Control 2	92 - 191	ng/mL
GAPDH	Control 1	1628 - 3381	pg/mL
	Control 2	6528 - 13557	pg/mL
Syn-1	Control 1	49 - 103	ng/mL
	Control 2	194 - 404	ng/mL
Flot-1	Control 1	32 - 65	ng/mL
	Control 2	148 - 307	ng/mL
TSG101	Control 1	2846 - 5910	ng/mL
	Control 2	11658 - 24212	ng/mL

Note: Argonaute-2 (Ago-2) is a Qualitative assay. No QC ranges are provided for this analyte.

HEXSM-6170: HEXSM-101.201 (15 December, 2023)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA