3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Vascular Endothelial Growth Factor, human recombinant, expressed in *E. coli*

Catalog Number **V7259** Storage Temperature –20 °C

CAS RN 127464-60-2

Synonyms: VEGF, vasculotropin

Product Description

Vascular Endothelial Growth Factor (VEGF) is an angiogenic growth factor, which is heat and acid stable. VEGFs stimulate endothelial cell growth, angiogenisis, and capillary permeability. VEGF is a secreted homodimeric, heparin-binding glycoprotein, which has an isoelectric point of 8.5. VEGF promotes the growth of endothelial cells isolated from bovine adrenal cortex, cerebral cortex, fetal and adult aorta, and human umbilical vein.

The target cell specificity of VEGF is restricted to vascular endothelial cells.² VEGF has no mitogenic effect on cultured corneal endothelial cells, vascular smooth muscle cells, BHK-12 fibroblasts, keratinocytes, human sarcoma cells, or lens epithelial cells.² A variety of human tumor cell lines including sarcoma and carcinoma cells show a 3.7 kb RNA transcript that hybridizes with the VEGF probe in a Northern blot.² Mouse sarcoma 180 cells express the VEGF mRNA and secrete a VEGF-like mitogen.³ Four cDNA clones, arising through alternative slicing and encoding mature human monomeric VEGF with 121, 165, 189, or 206 amino acids, have been identified.

This product is lyophilized from a sterile filtered buffered solution with no additives. It is a homodimer protein of the 165 amino acid splice variant, with a calculated molecular mass of 38.2 kDa.

Purity: ≥97% (SDS-PAGE)

The biological activity of human recombinant VEGF is measured by its mitogenic activity on human umbilical vein endothelial cells.

Endotoxin: ≤0.1 EU/μg protein

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Reconstitute the contents of the vial in water to a concentration of 0.1–1.0 mg/ml. After reconstitution, this solution can then be diluted into other aqueous buffers.

Storage/Stability

Store the product at -20 °C. The lyophilized protein is stable for up to a few weeks at 2–8 °C, but is best stored at -20 °C.

Reconstituted VEGF can be stored at 2–8 °C for up to one week. For prolonged use, freeze in working aliquots at –20 °C. Repeated freezing and thawing is not recommended.

References

- Ferrara, N., and Henzel, W.J., Pituitary follicular cells secrete a novel heparin-binding growth factor specific for vascular endothelial cells. Biochem. Biophys. Res. Commun., 161, 851 (1989).
- 2. Ferrara, N. et al., Molecular and biological properties of the vascular endothelial growth factor family of proteins., Endo. Rev., **13**, 18 (1992).
- 3. Clauss, M. et al., Vascular permeability factor: a tumor-derived polypeptide that induces endothelial cell and monocyte procoagulant activity, and promotes monocyte migration., J. Exp. Med., **172**, 1535 (1990).
- 4. Leung, D.W. et al., Vascular endothelial growth factor is a secreted angiogenic mitogen., Science, **246**, 1306 (1989).

DT,KAA,MAM 12/13-1