



MOUSE ANTI-INTEGRIN β 3 MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB2008
LOT NUMBER:	
QUANTITY:	100 μ g
SPECIFICITY:	Monoclonal antibody MAB2008 is specific for the β 3 subunit of the vitronectin receptor (VNR)J, which is expressed on endothelial cells, platelets, melanoma, and glioma cells. The antibody is of high affinity.
IMMUNOGEN:	Membranes of freshly isolated human platelets.
CLONE:	SAP
ISOTYPE:	IgG _{1k}
APPLICATIONS:	Immunohistochemistry (frozen and fixed tissue sections) Western blot Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	Human. Other species not tested.
FORMAT:	Purified immunoglobulin
PRESENTATION:	Liquid at 1 mg/mL in 0.02M phosphate buffer, 0.25M sodium chloride, pH 7.6, containing 0.1% sodium azide as a preservative.
STORAGE/HANDLING:	Maintain refrigerated at 2-8°C in undiluted aliquots for up to 12 months.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com