

Product Information

BAX, GST-tagged, human recombinant, expressed in *E. Coli* cells

Catalog Number **SRP5166**
Storage Temperature -70°C

Product Description

BAX is a proapoptotic protein of the BCL2 protein family. BAX forms a heterodimer with BCL2 and functions as an apoptotic activator. BAX interacts with and increases the opening of the mitochondrial voltage-dependent anion channel (VDAC), which leads to the loss in the mitochondrial membrane potential and the release of cytochrome c.¹ The expression of BAX gene is regulated by the tumor suppressor p53 and BAX has been shown to be involved in p53-mediated apoptosis. Multiple alternatively spliced transcript variants, which encode different isoforms, have been reported for this gene.²

Recombinant full-length human BAX was expressed in *E. coli* cells using an N-terminal GST tag. The gene accession number is NM_004324. Recombinant protein stored in 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 10 mM glutathione, 0.1 mM EDTA, 0.25 mM DTT, 0.1 mM PMSF, and 25% glycerol.

Molecular mass: ~49 kDa

Purity: 70–95% (SDS-PAGE, see Figure 1)


Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The product ships on dry ice and storage at -70°C is recommended. After opening, aliquot into smaller quantities and store at -70°C . Avoid repeated handling and multiple freeze/thaw cycles.

Figure 1.
SDS-PAGE Gel of Typical Lot
70–95% (densitometry)



References

1. Shimizu, S. et al., Bcl-2 family proteins regulate the release of apoptogenic cytochrome c by the mitochondrial channel VDAC. *Nature*, **399**, 483-487 (1999).
2. Meijerink, J.P.P. et al., Hematopoietic malignancies demonstrate loss-of-function mutations of BAX. *Blood*, **91**, 2991-2997 (1998).

DKF,MAM 10/11-1