




SIGMA QUALITY CONTROL TEST PROCEDURE

Product Information

Enzymatic Assay of β -N-ACETYLGLUCOSAMINIDASE (EC 3.2.1.30) from Jack Beans Sigma Prod. No. A-2264

PRINCIPLE:



Abbreviations:

PNP-NAG = p-Nitrophenyl N-Acetyl- β -D-Glucosaminide

NAG = N-Acetyl- β -D-Glucosamine

CONDITIONS: T = 25°C, pH = 5.0, $A_{400\text{nm}}$, Light path = 1 cm

METHOD: Spectrophotometric Stop Rate Determination

REAGENTS:

- A. 100 mM Citrate Buffer with 200 mM Sodium Chloride and 0.02% (w/v) Bovine Serum Albumin, pH 5.0 at 25°C
(Prepare 100 ml in deionized water using Citric Acid, Free Acid, Monohydrate, Sigma Prod. No. C-7129, Sodium Chloride, Sigma Prod. No. S-9625, and Albumin, Bovine, Sigma Prod. No. A-4503. Adjust to pH 5.0 at 25°C with 1 M NaOH.)
- B. 10 mM p-Nitrophenyl N-Acetyl- β -D-Glucosaminide Solution (PNP-NAG)
(Prepare 5 ml in deionized water using p-Nitrophenyl N-Acetyl- β -D-Glucosaminide, Sigma Prod. No. N-9376.)
- C. 200 mM Borate Buffer, pH 9.8 at 25°C
(Prepare 100 ml in deionized water using Boric Acid, Sigma Prod. No. B-0252. Adjust to pH 9.8 at 25°C with 1 M NaOH.)
- D. β -N-Acetylglucosaminidase Enzyme Solution
(Immediately before use, prepare a solution containing 0.05 - 0.1 unit/ml of β -N-Acetylglucosaminidase in cold Reagent A.)

**Enzymatic Assay of β -N-ACETYLGLUCOSAMINIDASE
(EC 3.2.1.30)
from Jack Beans
Sigma Prod. No. A-2264**

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable containers:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	0.40	0.40
Reagent B (PNP-NAG)	0.50	0.50

Mix by inversion and equilibrate to 25°C. Then add:

Reagent D (Enzyme Solution)	0.10	-----
-----------------------------	------	-------

Mix by inversion and incubate for exactly 10 minutes at 25°C. Then add:

Reagent C (Borate Buffer)	3.00	3.00
Reagent D (Enzyme Solution)	-----	0.10

Mix by inversion and transfer to suitable cuvettes. Record the $A_{400\text{nm}}$ for both the Test and Blank using a suitable spectrophotometer.

CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(\Delta A_{400\text{nm}} \text{ Test} - \Delta A_{400\text{nm}} \text{ Blank})(4)(\text{df})}{(10)(18)(0.1)}$$

4 = Total volume (in milliliters) of Assay

10 = Time of assay (in minutes) as per the Unit Definition

18 = Millimolar extinction coefficient of p-Nitrophenol at 400 nm¹

0.1 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

**Enzymatic Assay of β -N-ACETYLGLUCOSAMINIDASE
(EC 3.2.1.30)
from Jack Beans
Sigma Prod. No. A-2264**

UNIT DEFINITION:

One unit will hydrolyze 1.0 μ mole of p-nitrophenyl N-acetyl- β -D-glucosaminide to p-nitrophenol and N-acetyl- β -D-glucosamine per minute at pH 5.0 at 25°C.

FINAL ASSAY CONCENTRATION:

In a 1.00 ml reaction mix, the final concentrations are 50 mM citric acid, 100 mM sodium chloride, 0.01% (w/v) bovine serum albumin, 5 mM p-nitrophenyl N-acetyl- β -D-glucosaminide and 0.005 - 0.01 unit β -N-acetylglucosaminidase.

REFERENCES:

Bessey, O.A., Lowry, O.H., and Brock, M.J., (1946) *Journal of Biological Chemistry* **164**, 321-329

Li, S.-C. and Li, Y.-T. (1970) *Journal of Biological Chemistry* **245**, 5153-5160

Borooah, J., Leaback, D.H. and Walker, P.G. (1961) *Biochemical Journal* **78**, 106-110

NOTES:

1. The millimolar extinction coefficient is described in Bessey, O.A. et al. (1946).
2. This assay is based on the cited references.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.