

Product Information

E64 Ready Made Solution

Non-freezing solution, 1 mg/mL

SAE0154

Product Description

E-64 is an irreversible, potent and highly selective cysteine protease inhibitor. The trans-epoxysuccinyl group of E-64, its active moiety, irreversibly binds to an active thiol group of many cysteine proteases, such as papain, actinidase, and cathepsins B, H and L,^{2,3} to form a thioether linkage.

E-64 is reported to be one of the most effective low molecular weight inhibitors of trypsin-catalyzed hydrolysis.³ E-64 inhibited the activity of bleomycin hydrolase, and blocked the activity of a yeast cysteine protease gene (YCP1) which induces an increase in bleomycin metabolism.⁴ For *in vivo* studies, E-64 can be very useful, because it has specific inhibition, it is permeable in cells and tissues, it has low toxicity, it is easily synthesized, and it is stable.²

E-64 does not react with the functional thiol group of L-lactate dehydrogenase or creatine kinase, both of which are non-protease enzymes.^{5,6} E-64 does not inhibit serine proteases (except trypsin) like the cysteine protease inhibitors leupeptin and antipain.^{2,7} It does not react with low MW thiol compounds like 2-mercaptoethanol. E-64 has been used as an active site titrant.^{2,8,9} Synthetic¹ and natural⁵ methods for E-64 preparation have been reported.

Reagent

This E-64 Ready Made Solution is supplied as a proprietary 2.8 mM non-freezing solution formulation.

Storage/Stability

Storage at -20 °C is recommended. The product, as supplied, is stable for two years. For short time periods, the product can be stored at 2-8 °C.

Preparation Instructions

The stock solution may be diluted at a ratio of 1:100 to 1:1000 to achieve a working concentration in a range between 28 µM and 2.8 µM, respectively. The effective concentration range is 1-10 µM.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

References

1. Hanada, K. *et al.*, *Agric. Biol. Chem.*, **42**(3), 529-536 (1978).
2. Katunuma, N. and Kominami, E., *Methods Enzymol.*, **251**, 382-397 (1995).
3. Sreedharan, S.K. *et al.*, *Biochem. J.*, **316**(3), 777-786 (1996).
4. Pei, Z. *et al.*, *Mol. Pharmacol.*, **48**(4), 676-681 (1995).
5. Hanada, K. *et al.*, *Agric. Biol. Chem.*, **42**(3), 523-528 (1978).
6. Hashida, S. *et al.*, *J. Biochem.*, **88**(6), 1805-1811 (1980).
7. Barrett, A.J. *et al.*, *Biochem. J.*, **201**(1), 189-198 (1982).
8. Beynon, R.J. and Bond, J.S. (eds.), *Proteolytic Enzymes: A Practical Approach*. IRL Press (Oxford, UK), p. 244 (1989).
9. Buttle, D.J., and Barrett, A.J., *Biochem. J.*, **223**(1), 81-88 (1984).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

SAE0154pis Rev 02/22 EO,ARO,DT,SM,GCY

