

89166 Legionella Supplement (Legionella Growth Supplement; Legionella BCYE Growth Supplement)

An enrichment supplement recommended for enhanced growth of Legionella species

Composition
(per vial, sufficient for 500 ml medium)

Part A:
L-Cysteine hydrochloride 200.00 mg

Part B:
Ferric pyrophosphate, soluble 125.00 mg
Distilled water 5.00 ml

Directions

Rehydrate the contents of part A with 5 ml sterile distilled water. Warm up the refrigerated part B to room temperature and aseptically add Part A and Part B to 500 ml of sterile, molten, cooled (45-50°C) Buffered Charcoal Yeast Extract Agar, Base (Cat. No. 86558). Mix well and pour into sterile Petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

