

## shelf life information data sheet

### Product Family: Lynx® CDR

#### shelf life information

Expiration Date: 48 Months

Storage Conditions: Ambient Temperature

Applicable Catalog Number(s): CDRF1TN05, CDRF4HN05, CDRM3HN05, CDRM4HN05, CDRM6HN05, CDRM8HN05, CDRF3HN05, CDRF6HN05, CDRF8HN05

#### quality attributes studied

Bacterial Challenge

Leak Test

Burst Test

The stated expiration date designates the time prior to which a lot of the product is expected to remain within the approved shelf life specification if stored under defined conditions, and after which must not be put into use.

Not all products listed have undergone stability testing to determine the products' shelf lives. The design of the stability schedule is such that only select product(s) within a product family are tested. The design assumes that the products selected are representative of the family as they utilize the same materials of construction, production processes, and packaging components as the products within the family.

Updates are made to this data sheet on an as-needed basis; the date of signature on this document reflects the most recent revision.



24-APR-19

Thomas Brackett  
Supplier Quality Engineering Manager  
Third Party Finished Goods Quality

The life science business of Merck KGaA, Darmstadt, Germany  
operates as MilliporeSigma in the US and Canada.

The vibrant M, Millipore, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.  
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

