

79894 MKTTn Novobiocin Supplement (Muller-Kauffmann Tetrathionate Novobiocin Supplement; Novobiocin Selective Supplement)

Supplement recommended for selective enrichment and isolation of Salmonellae.

Composition (per vial sufficient for 1000 ml medium) Novobiocin 40.00mg

Directions

Rehydrate the content of 1 vial aseptically with 5 ml sterile distilled water and aseptically add to sterile, cooled (45-50°C Mueller Kauffmann Tetrathionat Novobiocin Broth Base (Cat. No. 89176)). Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

