3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

Monoclonal Anti-LRP1, Clone LRP1-11 produced in mouse, purified immunoglobulin

Product Number L2420

# **Product Description**

Monoclonal Anti-LRP1 (mouse IgG1 isotype) is derived from the hybridoma LRP1-11 produced by the fusion of mouse myeloma cells and splenocytes from BALB/c mice immunized with a synthetic peptide corresponding to a fragment of human LRP1 (GeneID: 4035), conjugated to KLH. The corresponding sequence is identical in mouse, rat, dog, and monkey. The isotype is determined by ELISA using Mouse Monoclonal Antibody Isotyping Reagents, Product Number ISO2. The antibody is purified from culture supernatant of hybridoma cells grown in a bioreactor.

Monoclonal Anti-LRP1 recognizes human, rat, and mouse LRP1. The antibody may be used in various immunochemical techniques including immunoblotting (~500 kDa) and immunofluorescence.

LDL receptor-related protein 1 (LRP1), also known as  $\alpha_2$ -macroglobulin receptor, is a large endocytic receptor widely expressed on surface and intracellular membranes of many cells. LRP1 is a member of the LDL receptor family, which is involved in several biological processes. LRP1 is synthesized in the endoplasmic reticulum as a single-chain 600 kDa type I transmembrane glycosylated precursor that is cleaved subsequently in the Golgi complex into two subunits of  $\sim$ 515 and 85 kDa. The 515 kDa N-terminal subunit remains attached to the membrane through noncovalent association with the smaller 85 kDa C-terminal subunit that contains the transmembrane and cytoplasmic domains.

LRP1 expression is essential for early embryonic development. LRP1 functions in lipoprotein metabolism, proteinase regulation, activation of lysosomal enzymes, cellular entry of viruses and toxins, cell signaling pathways, cell migration, and neurodevelopment. LRP1 recognizes at least 30 different ligands including lipoproteins, proteinases, proteinase-inhibitor complexes, ECM proteins, bacterial toxins, viruses, and various intracellular proteins.<sup>1-5</sup>

# Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide as a preservative.

Antibody concentration: ~1.0 mg/mL

#### **Precautions and Disclaimer**

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Storage/Stability

Store at –20 °C. For continuous use, store at 2–8 °C for up to one month. For extended storage, freeze at –20 °C in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discard if not used within 12 hours.

### **Product Profile**

Immunoblotting: a working antibody concentration of 1-2  $\mu$ g/mL is recommended using a whole extract of human U87 cells.

<u>Note</u>: In order to obtain best results in various techniques and preparations, it is recommended to determine optimal working dilutions by titration.

#### References

- Strickland, D.K. et al., J. Biol. Chem., 265, 17401-17404 (1990).
- 2. Herz, J. et al., EMBO J., 9, 1769-1776 (1990).
- 3. Herz, J., and Strickland, D.K., *J. Clin. Invest.*, **108**, 779-784 (2001).
- 4. May, P. et al., Ann. Med., 39, 219-228 (2007).
- 5. Lillis, A.P. et al., *Physiol. Rev.*, **88**, 887-918 (2008).

VS,ST,TD,KAA,PHC,MAM 04/19-1