

## 91603 Listeria mono Selective Supplement II

A selective supplement recommended by ISO Committee for the isolation of *Listeria* species.

### Composition

(per vial, sufficient for 500 ml medium)

Ceftazidime	10.0 mg
Nalidixic acid, sodium salt	10.0 mg
Amphotericin B	5.0 mg

### Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of 0.2 N Sodium hydroxide, further add 3 ml of sterile distilled water. Mix well and aseptically add it to 460 ml of sterile, molten L. mono Differential Agar Base (Cat. No. 77408) along with sterile contents of L. mono Enrichment supplement I (Cat. No. 03708) and sterile rehydrated contents of L. mono Selective Supplement I (Cat. No. 92301) or add in 470 ml of L. mono confirmatory Agar Base (Cat. No. 92302) along with sterile contents of L. mono Enrichment Supplement II (Cat. No. 15895) and sterile rehydrated contents of L. mono Selective Supplement I (Cat. No. 92301).

Mix well and pour into sterile petri plates.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

