



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Rat Thyroid Hormone Magnetic Bead Panel Catalog # RTHYMAG-30K

Control Catalog # LRTH-6030

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
T3	Control I	RTH-115	941 – 1954	pg/mL
	Control II	RTH-215	2012 – 4179	pg/mL
T4	Control I	RTH-115	15866 – 32953	pg/mL
	Control II	RTH-215	38516 – 79995	pg/mL
TSH	Control I	RTH-115	159 – 330	pg/mL
	Control II	RTH-215	485 – 1008	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

LRTH-6030: RTH-115.215 (23 March 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.