



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Human Sepsis Magnetic Bead Panel 1 Kit
Kit Catalog# HSP1MAG-63K
Control Catalog # LHSP-6063
Lot# SPS-114 and SPS-214

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
MIF	Control I	267 – 554	pg/mL
	Control II	2562 – 5321	pg/mL
ICAM-1	Control I	1481 – 3077	pg/mL
	Control II	17046 – 35404	pg/mL
Fas Ligand	Control I	268 – 556	pg/mL
	Control II	2903 – 6029	pg/mL
Soluble Fas	Control I	209 – 433	pg/mL
	Control II	2395 – 4975	pg/mL
VCAM-1	Control I	1268 – 2634	pg/mL
	Control II	13314 – 27651	pg/mL
Total PAI-1	Control I	256 – 531	pg/mL
	Control II	3044 – 6321	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

LHSP-6063: SPS-114.214 (09 September 2020)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com