

80548 TSC Agar Supplement

For the isolation and enumeration of vegetative forms as well as spores from *Clostridium perfringens* in food, clinical specimens and other material. *C. perfringens* colonies are fluorescent when illuminated with UV light (360 nm). Recommended by the "Schweizerisches Lebensmittelhandbuch", 2001, chapter 56.

Composition

Sucrose	5.0 g/l
D-Cycloserine	0.4 g/l
4-Methylumbelliferylphosphate disodium salt	0.1 g/l

Direction

Rehydrate 5.5 g in 50 ml distilled water. Filter sterilize the solution and aseptically add to 1 litre molten, sterilized TSC Agar (Prod No 93745; at approx. 50°C), mix.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

