



Quality Control Ranges
MILLIPLEX[®] MAP Human IGF Binding Protein Panel Magnetic Bead
Panel
Kit Catalog # HIGFBMAG-53K
Control Catalog# HIGFBP-6053

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
HIGFBP1	Control I	IGFBP-114	0.44 – 0.91	ng/mL
	Control II	IGFBP-214	2.0 – 4.1	ng/mL
HIGFBP2	Control I	IGFBP-114	4.2 – 8.8	ng/mL
	Control II	IGFBP-214	19 – 39	ng/mL
HIGFBP3	Control I	IGFBP-114	0.7 – 1.4	ng/mL
	Control II	IGFBP-214	3 – 7	ng/mL
HIGFBP4	Control I	IGFBP-114	9.1 – 19	ng/mL
	Control II	IGFBP-214	39 – 81	ng/mL
HIGFBP5	Control I	IGFBP-114	17 – 35	ng/mL
	Control II	IGFBP-214	69 – 144	ng/mL
HIGFBP6	Control I	IGFBP-114	2.1 – 4.4	ng/mL
	Control II	IGFBP-214	8.5 – 18	ng/mL
HIGFBP7	Control I	IGFBP-114	0.8 – 1.6	ng/mL
	Control II	IGFBP-214	3.2 – 6.7	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

HIGFBP-6053: IGFBP-114.214 (28 APRIL 2025)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.SigmaAldrich.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA