

1.10337

Supelco

sigmaaldrich.com


100 Tests

MQuant®

Peroxide Test Peroxid-Test



For IfU please scan code with
My M Safety App or go to
sigmaaldrich.com



-



+



0

100

200


400

600

800

1000

mg/l H₂O₂



1 s



30 s

If the alert zone (!)
turns blue in color see
Instruction for Use

Bei Blaufärbung
der Warnzone (!) s.
Gebrauchsanweisung



©2021

Store cold and dry • (2 - 8 °C) • Kalt und trocken lagern

Lot

min. shelf life (YYYY/MM/DD)