



MOUSE ANTI-CD44var (v6) MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB4073
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1.0 mg/mL
SPECIFICITY:	Specific for an epitope encoded by exon v6 on the variant portion of human CD44.
ISOTYPE:	IgG ₁
APPLICATIONS:	Immunohistochemistry on acetone fixed cryostat sections using the alkaline phosphatase-anti-alkaline-phosphatase (APAAP) or peroxidase anti-peroxidase (PAP) procedures. Antibody may be used at a concentration of 1 - 2 µg/mL. Immunohistochemistry on paraffin embedded sections after microwave treatment (1 x 5 min). Antibody may be used at a concentration of 1 - 2 µg/mL. Immunocytochemistry on cell smears. Western blot Flow cytometry (1:20) Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	Human
FORMAT:	Purified immunoglobulin.
PRESENTATION:	Liquid in PBS with no preservative.
STORAGE/HANDLING:	Maintain frozen at -20°C for 12 months from date of receipt. Avoid repeated freeze/thaw cycles.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.