

## Product Information

# Anti-Human Kappa Light Chain (Bound and Free) Antibody, Mouse Monoclonal

Clone KP-53, Purified from Hybridoma Cell Culture

**K4377**

## Product Description

Monoclonal Anti-Human Kappa Light Chain (mouse IgG1 isotype) is derived from the hybridoma produced by the fusion of mouse myeloma cells and splenocytes from an immunized mouse. Purified Bence Jones kappa light chain protein was used as the immunogen.

Monoclonal Anti-Human Kappa Light Chain is immuno-specific for kappa light chains (bound or free) as determined by ELISA. The antibody does not react with human lambda light chains.

## Reagents

Supplied as a solution in 0.01 M PBS, pH 7.4, containing 15 mM sodium azide as a preservative.

## Precautions and Disclaimer

Due to the sodium azide content, a Safety Data Sheet (SDS) for this product has been sent to the attention of the safety officer of your institution. Consult the SDS for information regarding hazards and safe handling practices.

## Storage/Stability

For continuous use, store at 2–8 °C. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

## Product Profile

A dilution of 1:10,000 was determined by indirect ELISA using kappa light chain at 10 mg/mL as the coating protein.

In order to obtain best results, it is recommended that each individual user determine their optimum working dilution by titration assay.

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2002–2025 Merck and/or its affiliates. All Rights Reserved.  
23285433 Rev11/25

