pesign for sustainability (pfs) scorecard



With our DfS scorecard, we drive sustainability improvement during the product development process through multiple product sustainability criteria divided into seven impact areas.

Cellvento® Modifeed Prime COMP



Chemically defined cell culture feed, highly concentrated, no pH adjustment in hydration, ambient storage of hydrated liquid possible.

| Impact areas | | Results |
|--------------|---------------------------|--|
| 6 | MATERIALS | No sodium hydroxide is required for pH adjustment, reducing materials needed and exposure to a strong base |
| 每 | SUPPLIERS & MANUFACTURING | 100% of key suppliers participate in our supplier sustainability program and more than 50% of product manufacturing costs are covered by suppliers with an Ecovadis rating of bronze or better. |
| (†) | PACKAGING | Packaging meets our standards for plastic sustainability and recyclability, except the foil-lined bag in the 50L pack size which is required for moisture protection but is not recyclable |
| S 3 | ENERGY & EMISSIONS | Storage at +4°C is not required for hydrated Modifeed Prime solutions, potentially saving up to 365 kWh per year. Longer hydration time for Modifeed Prime means an overhead stirrer would need to be used longer, but based on manufacturer data, would consume less that 1 kWh per 5L batch of media |
| | WATER | No change compared to baseline product in consideration of our DfS criteria |
| * | USABILITY & INNOVATION | No sodium hydroxide required for pH adjustment and hydrated media can be stored ambient for up to 30 days |
| | CIRCULAR ECONOMY | No change compared to baseline product in consideration of our DfS criteria |

Baseline product: Cellvento® 4Feed

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M and Cellvento are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.