



QUALITY CONTROL RANGES
MILLIPLEX® MAP Mouse Soluble Cytokine Receptor Magnetic Bead Panel
Catalog # MSCRMAG-42K
QC Catalog # MSCR6042 Lot # MSCR-108A and MSCR-208A

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
sCD30	Control 1	704 - 1462	pg/mL
	Control 2	6894 - 14318	pg/mL
sGP130	Control 1	668 - 1388	pg/mL
	Control 2	6579 - 13665	pg/mL
sIL-1RI	Control 1	351 - 728	pg/mL
	Control 2	3333 - 6922	pg/mL
sIL-1RII	Control 1	348 - 724	pg/mL
	Control 2	3448 - 7160	pg/mL
sIL-2Rα	Control 1	337 - 700	pg/mL
	Control 2	3376 - 7011	pg/mL
sIL-4R	Control 1	351 - 729	pg/mL
	Control 2	3374 - 7008	pg/mL
sIL-6R	Control 1	692 - 1436	pg/mL
	Control 2	6818 - 14160	pg/mL
sRAGE	Control 1	748 - 1554	pg/mL
	Control 2	7365 - 15297	pg/mL
sTNFR1	Control 1	136 - 283	pg/mL
	Control 2	1333 - 2768	pg/mL
sTNFR2	Control 1	364 - 755	pg/mL
	Control 2	3389 - 7038	pg/mL
sVEGFR1	Control 1	640 - 1329	pg/mL
	Control 2	6680 - 13873	pg/mL
sVEGFR2	Control 1	686 - 1426	pg/mL
	Control 2	6789 - 14099	pg/mL
sVEGFR3	Control 1	707 - 1469	pg/mL
	Control 1	6793 - 14108	pg/mL

MSCR6042: MSCR-108A.208A (10 August 2021)
EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.