



Click's Medium

Product Information

COMPONENT	C 5572 [1X] g/L	COMPONENT	C 5572 g/L
INORGANIC SALTS		VITAMINS	
CaCl ₂	0.1854	D-Pantothenic Acid•½Ca	0.002
MgCl ₂ •6H ₂ O	0.1649	Choline Chloride	0.002
KCl	0.4	Folic Acid	0.002
KH ₂ PO ₄	0.06	myo-Inositol	0.004
NaHCO ₃	1.35	Niacinamide	0.002
Na ₂ HPO ₄	0.0479	Pyridoxal•HCl	0.002
NaCl	8.0	Riboflavin	0.0002
AMINO ACIDS		Thiamine•HCl	0.002
Alanine	0.0356	OTHER	
Arginine•HCl	0.3175	Glucose	1.0
Asparagine•H ₂ O	0.06	Phenol Red•Na	0.01
Aspartic Acid	0.0532	Pyruvic Acid•Na	0.275
Cystine•2HCl	0.0782	ADD	
Glutamic Acid	0.0588	L-Glutamine	0.584
Glycine	0.03		
Histidine•HCl	0.1047		
Isoleucine	0.13		
Leucine	0.13		
Lysine•HCl	0.1812		
Methionine	0.0375		
Phenylalanine	0.08		
Proline	0.046		
Serine	0.042		
Threonine	0.12		
Valine	0.115		
Tryptophan	0.025		
Tyrosine•2Na	0.1189		