



VCX 500

500 Watts

Processing Volume: 200µl – 1L*



The **VCX 500** offers 500 watts of power to handle a range of liquid processing applications. This model shares the same microprocessor based operating system and programming screens as the VCX750 system.

This model offers the same programmability as larger **VCX** versions.

A standard ½" diameter probe with replaceable tip is included and many optional accessories are available to meet the needs of any application.

FEATURES:

- Programmable
- Variable power output
- 10 hour timer
- Pulse mode
- Wattage and Energy display
- Temperature monitoring
- Energy setpoint



(stand is optional)

Ordering information: Model no. **VCX 500**

500 watt ultrasonic processor with ½" diameter probe and tool kit.

SPECIFICATIONS

Power Output	500 watt maximum	Converter Cable	Part #201-0300
Frequency	20kHz	Length	6 ft. (1.8m)
Dimensions	H x W x D: 9.25 x 7.5 x 13.5" (235 x 190 x 340mm)		
Ultrasonic Converter	Part #CV334	Tool Kit	Part #381-0005
Diameter	2.5" (64mm)		Includes: 2 spanner wrenches and 1 open-end wrench
Length	7.25" (183mm)		
Standard Probe	Part #630-0220	Optional:	
Tip Diameter	0.5" (13mm) with replaceable tip	Stand with Clamp	Part #830-00459
Length	5.5" (139mm)	Sound Enclosure	Part #830-00427
Material	Ti6Al4V titanium alloy	Temperature Probe	Part #830-00060

*The actual processing volume of a specific probe or accessory is application specific.