

Product Information

Anti-EpCAM ICD antibody, Mouse monoclonal

Clone 4A7, purified from hybridoma cell culture

SAB4200473

Product Description

Anti-EpCAM ICD (mouse IgG1 isotype) is derived from the hybridoma 4A7 produced by the fusion of mouse myeloma cells and splenocytes from BALB/c mice immunized with a synthetic peptide corresponding to the intracellular domain of human EpCAM (GeneID: 4072), conjugated to KLH. The isotype is determined by ELISA using Mouse Monoclonal Antibody Isotyping Reagents (Cat. No. ISO2). The antibody is purified from culture supernatant of hybridoma cells grown in a bioreactor.

Anti-EpCAM ICD recognizes the intracellular domain (ICD) of human and mouse EpCAM. The antibody can be used in several immunochemical techniques including immunoblotting (~ 40 kDa), immunoprecipitation and immunofluorescence.

Epithelial cell adhesion molecule (EpCAM) is an epithelial transmembrane glycoprotein, composed of an ectodomain, a transmembrane domain and a cytoplasmic domain. EpCAM is frequently expressed in human malignancies, normal stem and progenitor cells, and at lower levels, in normal epithelia. In vivo and in vitro studies have shown that its intracellular domain (ICD) fragment translocates to the nucleus after cleavage from EpCAM ectodomain and associates with other molecules to bind DNA resulting in cell proliferation and tumor formation.¹⁻⁴

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide as a preservative.

Antibody Concentration: ~ 1.0 mg/mL

Precautions and Disclaimer

This product is for research use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze at -20 °C in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

Immunoblotting: a working concentration of 0.25-0.5 µg/mL is recommended using whole extracts of human SW480 cells.

Immunofluorescence: a working concentration of 5-10 µg/mL is recommended using human MCF7 cells.

Note: In order to obtain the best results using various techniques and preparations, we recommend determining optimal working dilutions by titration.

References

1. Maetzel, D., et al., *Nat. Cell Biol.*, **11**, 162-171 (2009).
2. Ralhan, R., et al., *BMC Cancer*, **10**, 331 (2010).
3. Carpenter, G., and Red Brewer, M., *Cancer Cell*, **15**, 165-166 (2009).
4. Went, P., et al., *Hum. Pathol.*, **35**, 122-128 (2004).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

SAB4200473dat Rev 06/21

For research use only. Not for use in diagnostic procedures.

The Merck logo, consisting of the word "MERCK" in a bold, red, sans-serif font.