



Eppendorf MiniSpin (MiniSpin +)

Z605220 and Z606235

* Rotor 5452 000 018 is included with the MiniSpin and 5453 000 011 with MiniSpin +

ml	Tubes per rotor	RCF	Rotor included with centrifuge Cat. No. / supplier ref.	Swing out or fixed	Suspension lid Cat. No. / model	Adapters Cat. No. / model	No. of adapters per rotor
1.5-2.0	12	12,100 (14,000)	Z605220 / 5452 000 018 (Z606235 / 5453 000 011)*	F	Z653411 / 5452 702 000	None	None
0.5-0.6	12	12,100 (14,000)	Z605220 / 5452 000 018 (Z606235 / 5453 000 011)*	F	Z653411 / 5452 702 000	Z366080 / 5425 716 001	12
0.4	12	12,100 (14,000)	Z605220 / 5452 000 018 (Z606235 / 5453 000 011)*	F	Z653411 / 5452 702 000	Z366099 / 5425 717 008	12
0.2	12	12,100 (14,000)	Z605220 / 5452 000 018 (Z606235 / 5453 000 011)*	F	Z653411 / 5452 702 000	A4190 / 5425 715 005	12