


pFLAG-CMV-5.1 (4.7 kb)



Multiple Cloning Site

(pFLAG-CMV-5.1)

Hind III	Not I	EcoR I	Cla I	Bgl II
ATT AAG CTT	GCG GCC	GCG AAT TCA	TCG ATA	GAT CTG
TAA TTC GAA	CGC CGG	CGC TTA AGT	AGC TAT	CTA GAC

EcoR V	Kpn I	Sal I	Xba I	BamH I	Sma I
ATA TCG GTA	CCA GTC	GAC TCT	AGA GGA	TCC CGG	GCT
TAT AGC CAT	GGT CAG	CTG AGA	TCT CCT	AGG GCC	CGA

FLAG Peptide Sequence

Asp	Tyr	Lys	Asp	Asp	Asp	Asp	Lys
-----	-----	-----	-----	-----	-----	-----	-----

GAT TAC AAG GAC GAC GAT GAC AAG TAG GTT

CTA ATG TTC CTG CTG CTA CTG TTC ATC CAA