



QUALITY CONTROL RANGES

MILLIPLEX® MAP Human Bone RANKL Single-Plex Panel Catalog # HRNKLMAG-51K-01

Control Catalog # LHBN6051RANKL

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
RANKL	Control I	RNKL-114	167 - 346	pg/mL
	Control II	RNKL-214	868 - 1803	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

LHBN6051RANKL: RNKL-114.214 (08 January 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA