

Technical Data Sheet

Anaerobic Jar

Ordering number: 1.13681.0001

Mode of Action

For cultivation of anaerobic and microaerophilic microorganisms in defined atmospheric conditions

- using Anaerocult® A (Cat. No. 1.32381.0001) and Anaerotest® (Cat. No. 1.32371.0001) for anaerobic microorganisms
- using Anaerocult® C (Cat. No. 1.32383.0001) for microaerophilic microorganisms.

The anaerobic jar is to be used together with Cat. No. 1.07040.0001 Petri dish rack (for up to 12 Petri dishes).

Ordering Information

Product	Cat. No.
Anaerobic Jar	1.13681.0001