

XXXXXX

7.76077.0007 01/2025

Supelco

sigmaaldrich.com

50 Tests

1.16997





Reflectoquant®

Total Hardness Test for RQflex® 20

Gesamthärte-Test für RQflex® 20

0.1 - 30.0 °d (0.1 - 38 °e 0.2 - 54 °f)

Reaction time: 15 sec



1°d \triangleq 10 mg/l CaO \triangleq 1.25 °e \triangleq 1.78 °f
 \triangleq 17.8 mg/l CaCO₃ \triangleq 0.18 mmol/l CaCO₃ (Ca)

Please consult Instruction for Use for further information and observe the RQflex® manual.



For IFU please scan code with My M Safety App or go to sigmaaldrich.com



©2025

Store cool and dry • (15 - 25 °C) • Kühl und trocken lagern

Lot

min. shelf life (YYYY/MM/DD)

Merck KGaA, 64271 Darmstadt, Germany
 EMD Millipore Corporation, Burlington, MA 01803, USA
 MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,
 Ontario, L6H 6J8, Canada