



QC RANGE SHEET
MILLIPLEX® MAP Human CVD Magnetic Bead Panel 6
Kit Catalog# HCVD6MAG-67K
QC Catalog# HCVMG6-6067-6
Lot# HCVM6-105 and HCVM6-205

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
CD14	Control 1	0.16 – 0.33	ng/mL
	Control 2	1.3 – 2.7	ng/mL
DPP4	Control 1	0.79 – 1.6	ng/mL
	Control 2	6.1 – 13	ng/mL
Endostatin	Control 1	0.04 – 0.08	ng/mL
	Control 2	0.30 – 0.63	ng/mL
Progranulin	Control 1	0.39 – 0.80	ng/mL
	Control 2	3.2 – 6.6	ng/mL
LRG1	Control 1	0.37 – 0.77	ng/mL
	Control 2	3.2 – 6.6	ng/mL
MCAM	Control 1	0.41 – 0.85	ng/mL
	Control 2	3.2 – 6.7	ng/mL
ucMGP	Control 1	0.59 – 1.3	ng/mL
	Control 2	5.3 – 11	ng/mL

HCVMG6-6067-6: HCVM6-105.205 (rev 17 July 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA