

Product Information

HumanKine™ Granulocyte-Colony Stimulating Factor, human recombinant, expressed in HEK 293 cells

Catalog Number **H5541**
Storage Temperature -20°C

Synonym: G-CSF

Product Description

HumanKine™ G-CSF is expressed as a glycosylated 21–25 kDa monomer in human 293 cells. Production in human 293 cells offers authentic glycosylation. Glycosylation contributes to stability in cell growth media and other applications.

Granulocyte CSF is produced by monocytes and fibroblasts.^{1,2} It stimulates granulocyte colony formation, activates neutrophils and mature granulocytes, and promotes differentiation of certain myeloid leukemic cells. Human and murine G-CSF have about 73% homology and show biological cross-reactivity.

Four distinct colony-stimulating factors (CSFs) that promote survival, proliferation, and differentiation of bone marrow precursor cells have been well characterized: granulocyte/macrophage-CSF (GM-CSF), granulocyte-CSF (G-CSF), macrophage-CSF (M-CSF), and interleukin-3 (IL-3, Multi-CSF).^{3,4} Both GM-CSF and IL-3 are multipotential growth factors, stimulating proliferation of progenitor cells from more than one hematopoietic lineage. In contrast, G-CSF and M-CSF are lineage-restricted hematopoietic growth factors, stimulating final mitotic divisions and the terminal cellular maturation of the partially differentiated hematopoietic progenitors.

This product is lyophilized from a PBS solution.

ED₅₀: ≤ 0.16 ng/mL

The specific activity was determined by the dose dependent stimulation of the proliferation of murine M-NFS-60 cells (Mouse Myeloid Leukemia indicator cell line).

Purity: $\geq 95\%$

Endotoxin level: ≤ 1 EU/ μg

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Briefly centrifuge the vial before opening. It is recommended to reconstitute the protein in sterile PBS containing 0.1% endotoxin-free recombinant human serum albumin.

Storage/Stability

Store the product at -20°C . The lyophilized product remains active for one year at -20°C .

Upon reconstitution, the cytokine can be stored at $2-8^{\circ}\text{C}$ for short term only, or at -20°C to -80°C in aliquots for long term. Avoid repeated freeze-thaw cycles.

References

1. Metcalf, D., *Cell*, **43**, 5 (1985).
2. Groopman, J., *Cell*, **50**, 5 (1987).
3. Mazur, E., and Cohen J., *Clin. Pharmacol. Ther.*, **46**, 250 (1989).
4. Morstyn, G., and Burgess, A., *Cancer Res.*, **48**, 5624 (1988).

HumanKine is a trademark of HumanZyme Inc.

GS,JF,MAM 12/10-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.