



QUALITY CONTROL RANGES
MILLIPLEX® MAP Mouse Soluble Cytokine Receptor Magnetic Bead Panel
Catalog # MSCRMAG-42K
QC Catalog # MSCR6042 Lot # MSCR-110 and MSCR-210

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
sCD30	Control 1	633 - 1315	pg/mL
	Control 2	6334 - 13156	pg/mL
sGP130	Control 1	615 - 1277	pg/mL
	Control 2	6248 - 12977	pg/mL
sIL-1RI	Control 1	333 - 692	pg/mL
	Control 2	3229 - 6706	pg/mL
sIL-1RII	Control 1	328 - 682	pg/mL
	Control 2	3261 - 6772	pg/mL
sIL-2R α	Control 1	315 - 655	pg/mL
	Control 2	3217 - 6681	pg/mL
sIL-4R	Control 1	336 - 698	pg/mL
	Control 2	3347 - 6951	pg/mL
sIL-6R	Control 1	647 - 1343	pg/mL
	Control 2	6462 - 13421	pg/mL
sRAGE	Control 1	606 - 1259	pg/mL
	Control 2	6079 - 12626	pg/mL
sTNFR1	Control 1	132 - 273	pg/mL
	Control 2	1323 - 2748	pg/mL
sTNFR2	Control 1	336 - 697	pg/mL
	Control 2	3361 - 6981	pg/mL
sVEGFR1	Control 1	648 - 1345	pg/mL
	Control 2	6599 - 13706	pg/mL
sVEGFR2	Control 1	631 - 1310	pg/mL
	Control 2	6344 - 13177	pg/mL
sVEGFR3	Control 1	611 - 1270	pg/mL
	Control 1	6248 - 12977	pg/mL

MSCR6042: MSCR-110.210 (12 February 2025)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA