



**CHICKEN ANTI-GOAT IGG (H & L)
PEROXIDASE (HRP) CONJUGATED
AFFINITY PURIFIED ANTIBODY**

CATALOG NUMBER:	AP163P	QUANTITY:	1 mg
LOT NUMBER:		CONCENTRATION:	1 mg/mL
		HOST/ISOTYPE:	Chicken IgG
SPECIFICITY:	Goat IgG (H&L)		
APPLICATIONS:	<u>Western Blot</u> : 1:1,000-1:10,000 colormetric 1:1,000-1:30,000 ECL <u>Immunohisto/cytochemistry</u> : 1:200-1:500 <u>ELISA</u> : 1:10,000-1:100,000 <i>Optimal working dilution must be determined by the end user.</i>		
SPECIES REACTIVITY:	Reacts with goat IgG. Reactivity with other species has not been determined.		
IMMUNOGEN:	Purified goat IgG.		
PRESENTATION:	Purified by affinity chromatography. Liquid in PBS containing 0.1% Pro-Clean 400.		
STORAGE/HANDLING:	Maintain refrigerated at 2°-8°C for up to twelve months from date of receipt. Avoid repeated freeze/thaw cycles.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2003-2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com