



QUALITY CONTROL RANGES

MILLIPLEX® MAP Human Sepsis Magnetic Bead Panel 1 Kit
Kit Catalog# HSP1MAG-63K
Control Catalog # LHSP-6063
Lot# SPS-117 and SPS-217

Standard Catalog # LHSP-8063
Lot# 8063-17K

Analyte		Standard Value
MIF		47,680 pg/mL

****Please see reverse side for Quality Control Ranges****

Analyte

QC Level

Expected Range

Units

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

LHSP-6063: SPS-117.217 (02 January 2025)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Human Sepsis Magnetic Bead Panel 1 Kit
Kit Catalog# HSP1MAG-63K
Control Catalog # LHSP-6063
Lot# SPS-117 and SPS-217

MIF	Control I	153 – 318	pg/mL
	Control II	1682 – 3494	pg/mL
ICAM-1	Control I	1581 – 3283	pg/mL
	Control II	17306 – 35944	pg/mL
Fas Ligand	Control I	293 – 608	pg/mL
	Control II	2904 – 6031	pg/mL
Soluble Fas	Control I	210 – 437	pg/mL
	Control II	2194 – 4557	pg/mL
VCAM-1	Control I	1253 – 2602	pg/mL
	Control II	12426 – 25808	pg/mL
Total PAI-1	Control I	252 – 523	pg/mL
	Control II	2948 – 6123	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

LHSP-6063: SPS-117.217 (02 January 2025)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA