



Specification

1.00012.0000 Acetone for gas chromatography ECD and FID SupraSolv®

Specification		
Purity (GC)	≥ 99.8	%
Identity (IR)	conforms	
Evaporation residue	≤ 3.0	mg/l
Water	≤ 0.05	%
Color	≤ 10	Hazen
GC/ECD (retention range 1,2,4-tri-chlorobenzene to decachlorobiphenyle, individual signals (lindane standard))	≤ 3	pg/ml
GC/FID (retention range n-undecane to n-tetracontane, individual signals (n-tetradecane standard))	≤ 3	ng/ml
suitable for residue analysis		

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.