



Specification

1.00264.0000 Formic acid 98-100% for analysis EMSURE® ACS, Reag. Ph Eur

Specification		
Assay (alkalimetric)	≥ 98.0	%
Color	≤ 10	Hazen
Acetic acid (CH ₃ COOH)	≤ 500	ppm
Chloride (Cl)	≤ 5	ppm
Sulfate (SO ₄)	≤ 5	ppm
Sulfite (SO ₃)	≤ 10	ppm
Heavy metals (as Pb)	≤ 10	ppm
Ag (Silver)	≤ 0.020	ppm
Al (Aluminium)	≤ 0.050	ppm
Ba (Barium)	≤ 0.050	ppm
Be (Beryllium)	≤ 0.020	ppm
Bi (Bismuth)	≤ 0.10	ppm
Ca (Calcium)	≤ 0.20	ppm
Cd (Cadmium)	≤ 0.050	ppm
Co (Cobalt)	≤ 0.020	ppm
Cr (Chromium)	≤ 0.050	ppm
Cu (Copper)	≤ 0.020	ppm
Fe (Iron)	≤ 2.0	ppm
Ge (Germanium)	≤ 0.050	ppm
K (Potassium)	≤ 0.10	ppm
Li (Lithium)	≤ 0.020	ppm
Mg (Magnesium)	≤ 0.50	ppm
Mn (Manganese)	≤ 0.050	ppm
Mo (Molybdenum)	≤ 0.020	ppm
NH ₄ (Ammonium)	≤ 10	ppm
Na (Sodium)	≤ 0.50	ppm
Ni (Nickel)	≤ 0.050	ppm
Pb (Lead)	≤ 0.020	ppm
Sr (Strontium)	≤ 0.020	ppm
Ti (Titanium)	≤ 0.10	ppm
Tl (Thallium)	≤ 0.050	ppm
V (Vanadium)	≤ 0.050	ppm
Zn (Zinc)	≤ 0.050	ppm
Zr (Zirconium)	≤ 0.10	ppm
Evaporation residue	≤ 10	ppm
Dilution test	passes test	
ACS, Ph Eur-Reagent		
Expiry date : see product label		

Specification

1.00264.0000 Formic acid 98-100% for analysis EMSURE® ACS,Reag. Ph Eur

Christian Herrmann

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.