



Specification

1.00397.0000 N,N-Dimethylformamide for peptide synthesis

Specification		
Purity (GC)	≥ 99.9	%
Identity (IR)	conforms	
Color	≤ 10	Hazen
Acidity	≤ 0.0005	meq/g
Alkalinity	≤ 0.0002	meq/g
Amines (Dimethylamine)	≤ 10	ppm
Evaporation residue	≤ 10	mg/l
Water	≤ 300	ppm

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.