



# Specification

---

1.01115.0000 Ammonium acetate EMPLURA®

---

| Spec. Values                     |             |   |
|----------------------------------|-------------|---|
| Assay (alkalimetric)             | ≥ 96.0      | % |
| Identity                         | passes test |   |
| pH-value (5 %; water)            | 6.0 - 7.5   |   |
| Chloride (Cl)                    | ≤ 0.002     | % |
| Sulfate (SO <sub>4</sub> )       | ≤ 0.01      | % |
| Heavy metals (as Pb)             | ≤ 0.0005    | % |
| Fe (Iron)                        | ≤ 0.001     | % |
| Residue on ignition (as sulfate) | ≤ 0.02      | % |
| Water                            | ≤ 2.5       | % |

Claudia Wiegand

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.