



# Specification

1.02378.0000 Calcium chloride anhydrous granules Reag. Ph Eur

Spec. Values		
Assay (calculated on dried substance)	≥ 98.0	%
Free alkali (as Ca(OH) <sub>2</sub> )	≤ 0.2	%
Water absorption capacity (24 hrs., 80 % relative humidity)	≥ 25.0	%
Loss on drying (200 °C)	≤ 5.0	%
Conforms Reag. Ph.Eur.		

Kristin Starke  
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.