



Specification

1.04169.0000 Glycine buffer substance for electrophoresis

Specification		
Assay (perchloric acid titration)	≥ 99.7	%
Identity (IR-spectrum)	passes test	
Chloride (Cl)	≤ 0.003	%
Heavy metals (as Pb)	≤ 0.001	%
NH ₄ (Ammonium)	≤ 0.02	%
Foreign amino acids	≤ 0.1	%
Other ninhydrine positive substances	≤ 0.1	%
Suitability (as buffer substance for electrophoresis)	passes test	

Dr. Timo Huxel

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.