



Specification

1.06443.0000 Sodium formate for analysis EMSURE® ACS, Reag. Ph Eur

Spec. Values		
Assay (perchloric acid titration)	≥ 99.0	%
Insoluble matter	≤ 0.005	%
pH-value (5 %; water)	7.0 - 8.5	
Chloride (Cl)	≤ 0.001	%
Phosphate (PO ₄)	≤ 0.001	%
Sulfate (SO ₄)	≤ 0.001	%
Ca (Calcium)	≤ 0.0005	%
Fe (Iron)	≤ 0.0005	%
Heavy metals (ACS)	≤ 0.0005	%
Loss on drying (150 °C)	≤ 0.5	%

Corresponds to ACS, Reag. Ph Eur

Kristin Starke

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.