



# Specification

1.07288.0000 n-Pentane hypergrade for organic trace analysis SupraSolv®

Specification		
Purity (GC)	≥ 99.9	%
Identity (IR)	conforms	
Evaporation residue	≤ 3.0	mg/l
Water	≤ 0.01	%
Color	≤ 10	Hazen
GC/ECD (retention range 1,2,4-tri-chlorobenzene to decachlorobiphenyle, individual signals (lindane standard))	≤ 2	pg/ml
GC/ECD (retention range dichloromethane to 1,2,4-trichlorobenzene individual signals (tetrachloromethane))	≤ 1	ng/ml
GC/FID (retention range n-undecane to n-tetracontane, individual signals (n-tetradecane standard))	≤ 2	ng/ml
GC/MSD (retention range n-undecane to n-tetracontane, scanning area 30 - 600 amu, individual signals (n- tetradecane standard))	≤ 2	ng/ml

suitable for residue analysis

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.