



Specification

1.13792.0001 PLC Silica gel 60 F₂₅₄, 1 mm with concentrating zone 20 x 4 cm 15 Glass plates 20 x 20 cm

	hRf-values
Colour test	
Bleu VIF Organol	17 - 27
Ceres black	41 - 51
Ceres violet BRN	57 - 72

Typical value determined on a conditioned plate
Eluent: Toluene (45% rel. humidity)

Tom Kupfer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.