



# Specification

1.50649.0001 Purospher® STAR RP-18 endcapped (2µm) Hibar® HR  
150-2.1 UHPLC column

	Spec. Values
Pressure	≤ 420 bar
Theoretical Plates (N/m) (1,2-Benzanthracene)	≥ 170000
Symmetry (1,2-Benzanthracene)	0.8 - 1.4
Capacity factor (1,2-Benzanthracene)	8.5 - 11.5

Peter Schaub

responsible laboratory manager quality control

*This document has been produced electronically and is valid without a signature*